

AED PADS

Compatible AED Defibrillation,
Stimulation and training

Pads for various manufacturers



**Adult Standard and
Radiolucent Pads**




Adult Pads Lead out




**Pediatric standard
Radiolucent Pads**

INTRO
Agilent, HP, Laerdal,
Philips


 = available

- EF 79XX : Adult Standard
- EF 79XX W: Adult Lead Out
- EF 77XX: Adult Radiolucent
- EF 77XX W: Adult Radiolucent + Lead Out
- EF 79XX P: Pediatric Standard
- EF 79XX PW: Pediatric Standard Lead Out
- EF 77XX P: Pediatric Radiolucent
- EF 77XX PW: Pediatric Radiolucent + Lead out
- EF 79XX PL: Pediatric Standard – Reduced Energy
- EF 79XXWPL: Pediatric Standard – Reduced Energy – Lead Out

Agilent – Hewlett Packard - Laerdal



Connectors	Model	EF7953	EF7953W	EF7753	EF7953P	EF7953PW	EF7753
	Codemaster 100						
	Codemaster XE						
	Codemaster XL						
	Codemaster XL+						
	Heartstream XL e- and XLT, with cable M3507A						
	Laerdal Heartstart 4000						

Agilent – Phillips (See also Phillips)


Connectors	Model	EF7950	EF7950W	EF7750	EF7750W	EF7950P (1)	EF7950PW (1)	EF7950PL	EF7950WP L	EF7750PL (1)
 <p>EF7950P, EF7750P, EF7950PW</p> <p>EF7950PL, EF7950WPL</p> <p>EF7950, EF7950W, EF7750</p>	Heartstream / Heartstart XL model 4735A with Philips M3508A Cable									
	Heartstream / Heartstart XLT model M3500B cable with Philips M3508A									
	Heartstream / Heartstart XL model 4735A with Philips M3508A Cable									
	Heartstart 4000 model 861290 with Philips cable M3508A									

EF7952PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J SEE NOTES !!!!




Artema / Innomed / S&W

Connectors	Model	EF7955	EF7955W	EF7755	EF7755W	EF7955P	EF7955PW	EF7755P
		 	ARTEMA CARDIOAID 360B ARTEMA CARDIOAID 200B ARTEMA CARDIOAID 200 ARTEMA CARDIOAID 100 ARTEMA CARDIOAID XL+ ARTEMA CARDIOAID XL S&W CARDIO AID LS S&W CARDIO AID MC S7W CARDIO AID MC+					

Cardiaid



Connectors	Model	EF7964W
	CARDIAID CT0207	

Cardioline


Connectors	Model	EF7952	EF7952W	EF7752	EF7752W	EF7952P (1)	EF7952PW (1)	EF7952PL	EF7952WPL	EF7752P (1)	EF7752PL (1)
		  EF7952 EF7952W EF7752 EF7752W EF7952P EF7952PW EF7752P EF7752PW  	ELIFE 700								

EF7952PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J SEE NOTES !!!!

Corpuls


Connectors	Model	EF7965/PLUS	EF7965	EF7965P	EF7965PXL
	Corpuls 3 NOT COMPATIBLE IF CORPULS 3 IS EQUIPPED WITH MODULAR SYSTEM DEFIB / PACER UNIT <i>SLIM</i>				

CU Medical

Connectors	Model	EF7958	EF79658W	EF7958P(1)	EF7958PL
	PARAMEDIC CU-ER1				
	PARAMEDIC CU-ER2 / ER-3				
	PARAMEDIC CU-ER5 WITH CABLE ER-5 OA3				


EF7958PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 150 J TO 50J SEE NOTES !!!

CU Medical


Connectors	Model	EF7959	EF7959W	EF7759	EF7959P(1)	EF7959PL	EF7959WPL
	I-PAD NF1200						

EF7959PL EF7959WPL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J see NOTES !!!


Defibtech

Connectors	Model	EF7966W	EF7966WPL
 EF7966W EF7966WPL	DDU – 100 Lifeline AED		
	DDU -120 Lifeline AUTO AED		


Drager

Connectors	Model	EF7955	EF7955W	EF7755	EF7755W	EF7955P	EF7955PW	EF7755P
	CARDIOLOG 2000							


GE

Connectors	Model	EF7981
	RESPONDER 3000 (DIRECT CONNECT)	
	RESPONDER 2000 (DIRECT CONNECT)	
	RESPONDER 1100 (DIRECT CONNECT)	
	CARDIOSERV (DIRECT CONNECT)	
	DASH RESPONDER (DIRECT CONNECT)	


Instramed

Connectors	Model	EF7984W	EF7984PW
	INSTRAMED AED		


Instramed

Connectors	Model	EF7985	EF7985 P
	INSTRAMED CARDIOMAX		


Innomed - STIM

Connectors	Model	EF7955S/ STIM
	CARDIAID 200B	


Marquette(GE)

Connectors	Model	EF7962	EF7762	EF7962P
	MARQUETTE 1250 FIRST RESPONDER			
	MARQUETTE 1500 FIRST RESPONDER			
	MARQUETTE RESPONDER 2000 (WITH ADAPTOR CABLE)			
	MARQUETTE RESPONDER 3000 (WITH ADAPTOR CABLE)			


Marquette (GE)

Connectors	Model	EF7972	EF7972 P
	CAVO – CABLE GE MEDICAL SYSTEMS REF: 22338301 SUPPLIED WITH THE FOLLOWING DEFIBRILLATORS – CARDIOSERV		





M & B (SEE ZOLL)

Connectors	Model	EF7951	EF7951W	EF7751	EF7751W
	AED 7000				

Mediana


Connectors	Model	EF7969	EF7969PW
 EF7969W EF7969PW	AED HeartOn A10		

Medtronic - Physiocontrol


Connectors	Model	EF7952	EF7952W	EF7752	EF7752W	EF7952P (1)	EF7952PW (1)	EF7952PL	EF7952WPL	EF7752P (1)	EF7752PL (1)
  EF7952 EF7952W EF7752 EF7752W EF7952P EF7952PW EF7752P EF7752PW   EF7952PL EF7952WPL	LIFEPAK 9 – connector QUICK COMBO										
	LIFEPAK 10										
	LIFEPAK 12										
	LIFEPAK 15										
	LIFEPAK 20										
	LIFEPAK 500 AED – WITH GREY CONNECTOR										
	LIFEPAK 500 AED – WITH PINK CONNECTOR										
	LIFEPAK 1000 AED										
	LIFEPAK CR PLUS AED										
	LIFEPAK EXPRESS AED										

EF7952PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J SEE NOTES !!!!


Medtronic - Physiocontrol

Connectors	Model	EF7971
	CAVO – CABLE PHYSIOCONTROL REF 3011030-01 (SUPPLIED IE: WITH THE FOLLOWING DEFIBRILLATORS – LIFEPAK 12)	


Medtronic - Physiocontrol

Connectors	Model	EF7973
	CAVO – CABLE PHYSIOCONTROL REF 802905-03	


Metrax Primedic

Connectors	Model	EF7960	EF7960W	EF7760
	HEARTSAVE AED			
	HEARTSAVE AED - M			
	HEARTSAVE 6			
	HEARTSAVE 6S			
	DEFIMONITOR XD			

Metsis

Connectors	Model	EF7983W	EF7983PW
	LIFEPOINT PRO		
	LIFEPOINT AED		

MRL

Connectors	Model	EF7955	EF7955W	EF7755	EF7755W	EF7955P	EF7955PW	EF7755P
	LITE							
	LIFEQUEST							

Mindray


Connectors	Model	EF7952	EF7952W	EF7752	EF7752W	EF7952P (1)	EF7952PW (1)	EF7952PL	EF7952WPL	EF7752P (1)	EF7752PL (1)
<p>EF7952 EF7952W EF7752 EF7752W EF7952P EF7952PW EF7752P EF7752PW</p> <p>EF7952PL EF7952WPL</p>	BENEHEART D3										
	BENEHEART D6										

EF7952PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J SEE NOTES !!!!


Nihon Kohden

Connectors	Model	EF7955	EF7955W	EF7755	EF7755W	EF7955P	EF7955PW	EF7755P
	CARDIOLIFE TEC-5500 SERIES WITH CABLE JC-755V							
	CARDIOLIFE TEC-7600 SERIES WITH CABLE JC-755V							
	CARDIOLIFE TEC-7700 SERIES WITH CABLE JC-755V							
	CARDIOLIFE TEC-2202 WITH CABLE JC-755V							
	CARDIOLIFE TEC-2203 WITH CABLE JC-755V							
	CARDIOLIFE TEC-7531 WITH CABLE JC-755V							
	CARDIOLIFE TEC-8251 WITH CABLE JC-755V							
	CARDIOLIFE TEC-8300 SERIES WITH CABLE K342A(JC855V)							


Nihon Kohden

Connectors	Model	EF7961	EF7961W	EF7761	EF7961P
	CARDIOLIFE TEC-2200 SERIES – WITH CABLE JC-765V				
	CARDIOLIFE TEC-2300 SERIES – WITH CABLE JC-765V				
	CARDIOLIFE TEC-5500 SERIES – WITH CABLE JC-765V				
	CARDIOLIFE TEC-7600 SERIES – WITH CABLE JC-765V				
	CARDIOLIFE TEC-7700 SERIES – WITH CABLE JC-765V				
	CARDIOLIFE TEC-7531 – WITH CABLE JC-765V				
	CARDIOLIFE TEC-8300 SERIES – WITH CABLE K342B (JC-865V)				



Nihon Kohden

Connectors	Model	EF7963
	TEC-8250	
	TEC-8251	

Nihon Kohden - STIM




Connectors	Model	EF7955 / STIM
	TEC-7731 WITH CABLE JC-755V	
	TEC-5531 WITH CABLE JC-755V	

Osatu Bexen

Connectors	Model	EF7952	EF7952W	EF7752	EF7752W	EF7952P (1)	EF7952PW (1)	EF7952PL	EF7952WPL	EF7752P (1)	EF7752PL (1)
		 EF7952 EF7952W EF7752 EF7752W EF7952P EF7952PW EF7752P EF7752PW  EF7952PL EF7952WPL	REANIBEX – 200								
	REANIBEX – 300										
	REANIBEX – 500										
	REANIBEX – 700										
	REANIBEX - 800										

EF7952PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 200 J TO 50J SEE NOTES !!!!


Phillips / Laerdal (See also AGILENT)

Connectors	Model	EF7950	EF7950W	EF7750	EF7750W	EF7950P (1)	EF7950PW (1)	EF7950PL	EF7950WPL	EF7750PL (1)
		 EF7950P, EF7750P, EF7950PW  EF7950PL, EF7950WPL  EF7950 EF7950W, EF7750	Heartstream ForeRunner I AED model E, S, EM							
	Heartstart FR1 AED models 940010xx, 940020xx, 940030XX									
	Heartstream Heartstart FR2 AED models M3841A M3861A									
	Heartstream Heartstart FR2 + AED with ECG module model M3840A									
	Codemaster models – M1722A/B , M1723A/B, M1724A, M2475B – with cable Philips M3507A and pad adaptor 05-200									
	Heartstart MRx model M3535A – with cable Philips M3508A)									


EF7950PL REDUCED ENERGY : DEFIBRILLATOR OUTPUT IS REDUCED FROM 150J TO 50J

See NOTES !!!


Schiller

Connectors	Model	EF7954	EF7954W	EF7754	EF7954P	EF7754P
	FRED					
	FRED TELE					
	DEFIGARD DG 1002					
	DEFIGARD DG 3002					
	DEFIGARD DG 2000 EDOS					


Schiller

Connectors	Model	EF7956	EF7956W	EF7756	EF7956P
 <p>EF7956 EF7956W EF7956P EF7756</p>	FRED				
	FRED EASY 2G				
	DEFIGARD DG 4000				
	DEFIGARD DG 5000				
	DEFIGARDS DG 6002				

Schiller


Connectors	Model	EF7967	EF7967P
 <p>EF7967 EF7967P</p>	FRED EASYPORT		
	ARGUS PRO Lifecare		
	ARGUS PRO Lifecare 2		
	CARDIOVIT AT101 Easy		

Welch Allyn


Connectors	Model	EF7955	EF7955 W	EF7755	EF7755 W	EF7955 P	EF7955 PW	EF7755 P
	AED - 10 WITH PAEDIATRIC ADAPTOR WA 002173-U (2)							
	AED - 20 WITH PAEDIATRIC ADAPTOR WA 002173-U (2)							
	PIC 30							
	PIC 40							
	PIC 50							

(2) IN PEDIATRIC APPLICATIONS – THE USE OF PLATES WITHOUT THE INDICATED ENERGY REDUCER ADAPTOR CAN BE DANGEROUS


Weinmann

Connectors	Model	EF7964 W
	MEDUCORE EASY	





Zoll




Connectors	Model	EF7951	EF7951W	EF7751	EF7751W	EF7951P	EF7951PW	EF7751P
	PD 1200							
	PD 1400							
	PD 1600							
	PD 1700							
	PD 2000							
	M – SERIES							
	E – SERIES							
	R - SERIES							

Zoll

Connectors	Model	EF7951 / STIM	EF7951P / STIM
	M – SERIES (PACER VERSIONS ONLY)		

ADAPTOR LISTINGS

Connectors	MANUFACTURES	EFXX 53
	AGILENT HEWLETT PACKARD LAERDAL	
Connectors	MANUFACTURERS	EFXX 50
 <p>EF7950P, EF7750P, EF7950PW</p>  <p>EF7950PL, EF7950WPL</p>  <p>EF7950 EF7950W, EF7750</p>	PHILIPS AGILENT LAERDAL	

Connectors	MANUFACTURERS EFXX 55
	<p>ARTEMA INNOMED S&W Drager MRL Nihon Kohden Welch Allyn</p>
Connectors	MANUFACTURERS EFXX 52
 <p>EF7952 EF7952W EF7752 EF7752W EF7952P EF7952PW EF7752P EF7752PW</p> <p>EF7952PL EF7952WPL</p>	<p>Medtronic Physiocontrol Osatu Bexen Cardioline Mindray</p>
Connectors	MANUFACTURERS EFXX 64
	<p>Cardaid Weinmann</p>

NOTES

- (1) Do not use with fully automatic defibrillators (AED) ; in semi- automatic defibrillators use only in manual mode
- (2) In pediatric applications – the use of plates without the indicated energy reducer adaptor can be dangerous