

EvoCare[®]

ECG

Compatible patient
cables with 10
leads, banana,
needle, snap or clip

One-Piece Cables

3.6m to limbs
3.2m to chest



Leadwire Sets

1.2m to limbs 0.8m to
chest



Trunk Cables


2.4m cable




ECG Accessories



Choosing your ECG patient cable




10 Leads Banana (4mm)
Used together with limb clamps or Suction electrodes.
With Adaptors for electrodes or snap electrodes



10 Leads Needle (3mm)
NOT AVAILABLE



10 Lead Snap
Moulded Snap termination
Are used with press stud electrodes. Select models only



10 Lead Clip
Moulded Grabber termination
are used with press stud electrodes. Select models only

Biocare One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
ECG 101/ 101G/ 300 / 300G/ 1200 9803		10 Leads Banana	EV10-BK1-B/I	012-0700-00, 007517, 007514, 007705, 007706	EVA002BNI
		10 Leads Clip	EV10-BK1-P/I		EVA002CNI

Burdick One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
EK10, Elite, E350, E350I, E550, E560, Eclipse LE,		10 Leads Banana	EV10-BK1-B/I	012-0700-00, 007517, 007514, 007705, 007706	EVA002BNI
		10 Leads Clip	EV10-BK1-P/I		EVA002CNI
Eclipse 4, 300, 400, 800, 850, Plus and Premier, Elite II, Eclipse LE II		10 Leads Banana	EV10-BK2-B/I	007704, 007725, 012-0844-00, 012-0844-01, 007785	EVA004BNI
		10 Leads Clip	EV10-BK2-P/I		EVA004CNI

Dong Jiang One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
ECG 11A / 11B/ 11D/ 32A/ 1220		10 Leads Banana	EV10-BK1-B/I	012-0700-00, 007517, 007514, 007705, 007706	EVA002BNI
		10 Leads Clip	EV10-BK1-P/I		EVA002CNI





Cardioline (See Nihon Kohden)

Carewell (See Nihon Kohden)



Edan One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Edan SE-3, SE-12, SE-12 Express, SE-601A, SE-601B, SE-601C, SE-1010PC, SE-1200, SE-12 Express, ST-1212 Nihon Kohden BSM-2301K Spacelabs SL6, Mindray DECG-03A,		10 Leads Banana	EV10-ED-B/I	01.57.107048	EVA030BBI
		10 Leads Clip	EV10-ED-P/I		EVA030CBI





Fakuda Denshi One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
FX-101, FX-102F, FX-3010, FX-328UY, FCP 2155, FCP 2201, FCP 4101, FCP 4101A, CardiMax FX-7101, FX-7102, FX-7202, FX-7402	 	10 Leads Banana	EV10-FD1-B/I	CP-104L	EVA017BAI
		10 Leads Clip	EV10-FD1-P/I		
	 	12 Leads	EV-21099S		EVA014SNI
		12 Leads Clip	EV-21099P		EVA014CNI





Fakuda ME One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
KP-500, KP-500D, Cardisuny	 	10 Leads Banana	K10-FD2-B0		EVA024BNA
		10 Leads Clip	K10-FD2-P0		EVA024CEI 20KoHM+ Discharge
		10 Lead Snap			EVA024SEI 20KOHM +discharge

Hellige One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
EK36, EK41, EK43, EK53, EK56, EK403 EK413, EK512, CARDIOSMART, CARDIOSYS, CORONA	 	10 Leads Banana	EV10-HL-B/I		EVA029BNI
		10 Leads Clip	EV10-HL-P/I		EVA029CNI
EK51	 	10 Lead Banana			EVA011DBI

Kenz One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
103, 106, 107, 302, 303, 1201 and Cardioline Plus, ECG 106	 	10 Leads Banana	EV10-KZ1-B/I	PC-104, 63050075, 63050074, 895.0586	EVA016BNI
		10 Leads Clip	EV10-KZ1-P/I		EVA016CNI
ECG 108, 110, Cardico 302, 601, 1210, 1211	 	10 Leads Banana	EV10-KZ2-B/I		EVA015BNI
		10 Leads Clip	EV10-KZ2-P/I		EVA015CNI

Marquette (GE) One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
MAC 400, MAC 500, MAC 1000, MAC 1100, MAC 1200, MAC 1600, Microsmart M		10 Leads Banana	EV10-MQ-B/I		EVA008BBI
		10 Leads Clip	EV10-MQ-P/I		EVA008CBI
Eagle, Solar, Dash, MAC Lab, Cath Lab and Datex Ohmeda E Prestin module		10 Leads Snap	EV21086S/I		EV21086S/I
		10 Leads Clip	EV21086P/I		EV21086P/I

Marquette (GE) ECG Trunk Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
MAC 400, MAC 500, MAC 1000, MAC 1100, MAC 1200		10 Leads IEC	EV10-MQ	22341809	EVE008BAI
Eagle, Solar, Dash, Mac Lab, Cath Lab, and Datex Ohmeda E-Prestin and M-Prestin Modules		10 Leads IEC	EVQ-21086	2017006-001, 416035-001, 416035-003	EVE006BAI

Marquette (GE) ECG Leads

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
AM4 Acquisition Module on monitors including MAC PC, MACVU, MAC 6, MAC 8, MAC 12, MAC 15, CASE 12, CASE 15, CASE 16, MAX-Personal, MAX-1, CENTRA Systems with AM4 or AM5 Acquisition Modules		10 Leads Banana	EQD10-LB-I	900177-203	
		10 Leads Clip	EQD10-LP-I		EVF008CNI
		10 Leads Snap	EQD10-LS-I		EVF008SNI
MAC 400, MAC 500, MAC 1000, MAC 1100, MAC 1200		10 Leads Banana	EQ10-LB-I	38401816, 38401817	
		10 Leads Snap	EQ10-LS-I		EE051S5I
		10 Leads Clip	EQ10-LP-I		EE051C5I

Marquette (GE) ECG Leads

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
CAM 14 Acquisition Module or with a cable that requires CAM 14 Leadwires such as GE MAC 1100, 1200, 1200ST, MAC 1600, MAC 3500, MAC 5000, MAC 5000ST, MAC 5500, CASE Series, CASE 8000, CardioSoft, CardioSmart, CardioSys. Please note that these leadwires come with the option of no banana connectors (for use with the GE CAM 14 Electrode Connector Set 900178-001) or molded-in banana, snap, clip or needle connectors		10 Leads Banana	EQB10-LB-I		VE008BNI
		10 Leads Clip	EQB10-LP-I		VE008CNI
		10 Leads Snap	EQB10-LS-I		VE008SNI

Medtronic One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20		10 Leads Snap	EV21027S-I		VA013SBI
		10 Leads Clip	EV21027P-I		EV21027P-I



Mindray One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Mindray DECG-03A, Edan SE-3, SE-12, SE-12 Express, SE-601A, SE-601B, SE-601C, SE-1010PC, SE-1200, SE-12 Express, ST-1212 and Nihon Kohden BSM-2301K		10 Leads Banana	EV10-ED-B	01.57.107048	EV10-ED-B
		10 Leads Clip	EV10-ED-P		EV10-ED-P
Beneview 12Lead module		10 lead Snap		0010-30-42722	VA020SBI
		10 Lead Clip			VA020CBI



Mindray Trunk cable

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Mindray DECG-03A, Edan SE-3, SE-12, SE-12 Express, SE-601A, SE-601B, SE-601C, SE-1010PC, SE-1200, SE-12 Express, ST-1212 and Nihon Kohden BSM-2301K			EV10-ED-B	01.57.107048	EV10-ED-B


Mindray Trunk cable (Continue)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Beneview 12Lead module		10 lead IEC		01.57.107048	EVE027-BAI
Beneview 12Lead module MULTILINK		10 Lead			EVR013-BAI Used with VR024S6I AND VB051S4I


Mindray ECG Leads

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
All Multilink plugs		10 Leads Banana	EQB10-LB-I		VE008BNI
		10 Leads Clip	EQB10-LP-I		VE008CNI
		10 Leads Snap	EQB10-LS-I		VE008SNI
		10 Leads Snap			VR024S6I + VB051S4I



Mortara One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
ELI 10, ELI 50, ELI 100, ELI 150, ELI 150C, ELI 200, ELI 210, ELI 230, ELI 250, ELI 250C and ELI 350		10 Leads Banana	EV10-MT-B/I	9293-032-50, 9293-040-50	EV10-MT-B/I
		10 Leads Clip	EV10-MT-P/I	9293-042-50	EV10-MT-P/I
		10 Leads Snap	EV10-MT-S/I	9293-033-50	EV10-MT-S/I



Nihon Kohden One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Cardiofax, 6151, 6851, 6353, 8110, 8330, 8330A, Escape 12, Edan Smart ECG, Danatech, and Dr. Lee EKGs 120B		10 Leads Banana	10-NK1-B/I		VA001BAI
		10 Leads Clip	10-NK1-P/I		VA001CAI


Nihon Kohden One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Cardiofax Q, 1350, ECG-1250A, ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, ECG-9110K, ECG-9130K, Cardioline (new line), Del Mar Reynolds Voyager 12; Newtech 1206; Carewell 1103, Carewell AR 1200ADV		10 Leads Banana	EV10-NK2-B/I	BA-902D	VA009BNI
		10 Leads Clip	EV10-NK2-P/I	BA-903D	VA008CNI
		10 Lead Snap	EV21009S-I		VA036SNI
		10 Lead Clip	EV21009P-I	BJ-900P	VA036CNI

Philips HP One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
M1770A, M1771A, M1772A, PageWriter 100, 200, 200i and Agilent 300i, 300PI; CardioTech CRT-6U, CRT-6D, CRT-6T; Bionet Cardiocare 2000; Delmar Reynolds CardioCollect 12; BLT-08 LT ECG; Spacelabs CardioDirect 12; Woodway Cortex CPX System; Welch Allyn CP 100 and CP 200		10 Leads Banana	EV10-HP-B	M2461A, M3702C, M3703C, 400293 + 401129, 45502-13	VA018BCI
		10 Leads Clip	EV10-HP-P		VA018CNI
Codemaster 100		10 Leads Snap	EVP-21085S		P-21085S
		10 Leads Clip	EVP-21085P		P-21085P

Philips HP ECG Trunk Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
78352C, 78354C, 78834C Legacy HP Monitors, Modules C12 or Opt. C18, IntelliVue Multi-Modules M3002A, Opt. C12, Opt. C12, MP2 Monitor, M8105A Monitor, M3535A/M3536A, Defibrillator/Monitor. For use (LPA5-90S0) and M1602A (LPB5-90S0) lead set		10 Leads IEC	EVP-21085	M1663A, M1949A, M3525A,	EVP-21085

Philips HP ECG Leads

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Trim, Trim I, Trim II (989803130101), Trim III, XLi M1700A, PageWriter XLE Machines		10 Leads Banana	EVW10-LB	989803129161, M1713B, M1716B, M1702A	VX018BNI
		10 Leads Clip	EVW10-LP		EVW10-LP
		10 Leads Snap	EVW10-LS		EVW10-LS

Physio Control (See Medtronic)

Schiller One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
AT-1, AT-2, AT-2 Light, AT-2 Plus, AT-4, AT-10, AT-10 Plus, AT-60, AT-101, AT-102, AT-102-Plus, MS-2015; Bionet: CardioCare Series 40893, Biosound MyLab 30 Cardiette: AR1200 adv; Welch Allyn CP10, CP20		10 Leads Banana	EV10-SH1-B	2.400095, 2.400071S, 2.40071E, MD071, 80130-0000, 9293-032-60	VA022BNI
10 Leads Clip		EV10-SH1-P		VA022CNI	
AT3, AT6, CS6, AT5, AT10, AT60, CS-200; Bionet Cardiotouch		10 Leads Banana	EV10-SH2-B		VA021BNI
		10 Leads Clip	EV10-SH2-P		VA021CNI

Welch Allyn ECG Leads (Also see Schiller)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with the 400293 10-lead patient cable on models CP100, CP200		10 Leads Banana	EVWA10-LB		EVWA10-LB
		10 Leads Snap	EVWA10-LS		EVWA10-LS
		10 Leads Clip	EVWA10-LP		EVWA10-LP
Welch Allyn CP10, CP20		10 Leads Banana	EV10-SH1-B		VA021BNI
		10 Leads Snap	EV10-SH1-N		
		10 Leads Clip	EV10-SH1-P		VA021CNI

Zoll One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
1 Step Patient Cable for E Series and M Series defibrillators. Includes trunk cable and limb and V leads		10 Leads Snap	EV21069S	8000-1007-01, 8000-1007-02	VAO34SBI
		10 Leads Clip	EV21069P		VAO34CBI
		10 Leads Banana			VAO34BBI

ELECTRODES CLIPS AND ACCESSORIES

Connector	Options	EVOTEC OEM
	<p>SET OF 10 UNIVERSAL CLIPS</p>	<p>C01</p>
	<p>SET OF 10 GRABBERS TO 4MM BANANA</p>	<p>C02</p>
	<p>SET OF 10 GRABBERS WITH LEADS TO 4MM BANANA</p>	<p>V0010D</p>
	<p>SET OF 10 MULTI FUNCTION TO 4MM BANANA</p>	<p>C03</p>
	<p>SET OF 10 MULTI FUNCTION TO 4MM BANANA</p>	<p>C04</p>
	<p>CABLE SUPPORT TO KEEP CABLES OF FLOORS</p>	<p>C10</p>
	<p>UNIVERSAL ECG LIMB ELECTRODES IEC</p>	<p>C15</p>
	<p>ELECTRODE BULBS SET OF 6 TO 4MM BANANA</p>	<p>C16</p>