



EvoCare®

NIBP

Hoses,
reusable and
disposable
cuffs and
connectors

NIBP Hoses

8 foot / 2.4m length



NIBP Cuffs

All patient sizes, reusable and
disposable



NIBP Connectors





Left Blank

Bionet (See Quick-Release NIBP Hose)

Criticare (See Nihon Kohden and Quick-Release NIBP Hose)

Colin NIBP Hose (Also see Nihon Kohden)



Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
		Adult NIBP Hose	ES-31-32	A015ZZ	xxx
		Neonate NIBP Hose	ESN-31-20		

Corometrics (See Marquette)



Critikon (See Nihon Kohden)

Datascope (See Nihon Kohden and Quick-Release NIBP Hose)



Dinamap NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Dinamap ProCare, Pro Series, Pro 1000 and Critikon Compact BP		Adult NIBP Hose	ED-25-09	107365	Ya24A10-12
		Neonate NIBP Hose	EDN-25-27	107368	

Datex Ohmeda NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Datex Ohmeda AS/3, CS/3, S/5, Cardiocap 5, Cardiocap II, Light		Adult NIBP Hose	ED-22-17/18	877235	YA22A12-12
		Neonate NIBP Hose	EDN-22-27	877514	YA22A12-12S


Drager NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Siemens Drager SC6000, SC7000, SC8000, SC9000, Infinity Delta, Vista, Kappa, Gamma		Adult NIBP Hose	ES-40-15	MP00953, 1275 275 EH5OU, 1275 275 E07	XXX
		Neonate NIBP Hose	ESN-40-20	2870298	XXX


Edan (See Quick-Release NIBP Hose)

Goldway (See Quick-Release NIBP Hose)

Iris NIBP Hose



Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Omni 2 and Omni 3 Creative, Neosoft		Adult NIBP Hose	ES-36-2T		YA18A01-10

Marquette (GE) NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Dash, Eagle, Solar, Tram 400, Tram 450, Tram 451, Tram 800, Tram 850, Tram 851, Corometrics 120 and Corometrics 250, Carescape B850		Adult NIBP Hose	ED-24-17	414873-001, 9461-217, 2020980-001, 2017008-001	YA23A12-12
		Neonate NIBP Hose	EDN-24-27	414874-001	YA23A16-08
Eagle 4000, Solar, Tram 100, Tram 200 and Tram 300		Adult NIBP Hose	ED-22-17	9461-203	YA22A12-12
		Neonate NIBP Hose	EDN-22-27	877514	YA22A12-12N
Eagle 3000, Dash 2000, Dash 3000, Dash 4000, Dash 5000 and Corometrics 120		Adult NIBP hose	Xxx	Xxx	YA2310-12 A01
Eagle 3000, Dash 2000, Dash 3000, Dash 4000, Dash 5000 and Corometrics 120		Adult NIBP Hose	ED-24-09	2017008-003	YA2310-12
		Neonate NIBP Hose	EDN-24-20		YA23A16-08

**Medtronic, MEK,
Mennen, Mindray,
Nellcor, Nihon Kohden**

Medtronic NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Lifepak 12 and Lifepak 20		Adult NIBP Hose	ES-12-15		Y23A13-15
		Neonate NIBP Hose	ESN-12-20		YA23A13-15N




MEK (See Welch Allyn)

Mennen (See Welch Allyn)

Mindray (See Quick-Release NIBP Hose)

Nellcor (See Quick-Release NIBP Hose)

Nihon Kohden NIBP Hose

Models	Hose	Options	Unimed/Evo	ORIGINAL OEM #	EVOTEC OEM
		Adult NIBP Hose	ED-09-09	8841	YA14A10-12
		Neonate NIBP Hose	EDN-09-27	8840	YA14A10-12N
Life Scope OPV 1500K, BSM-2351A, BSM-4102, BSM- Scope TR, BSM-6000 series, BSM- BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701		Adult NIBP Hose	ED-30-18	YN-900P, YN-901P	YA28A11-08
		Neonate NIBP Hose	END-30-20	YN-920P, YN-921P	YA28A11-08N

**Nihon Kohden, Philips,
 Physio Control, Quick
 Release.**

Nihon Kohden NIBP Hose (Continued)

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Nihon Kohden, Criticare: 507E/P/R, 507N/P, 507S/SD, 5070/P, 1100, 2200; Datascope, Colin, Welch Allyn		Adult NIBP Hose	ES-04-07	709	YA18A01-10
		Neonate NIBP Hose	ESN-04-20		YA18A01-10N

Philips NIBP Hose (Also see Quick-Release NIBP Hose)

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
M1008B and all Philips monitors except for A1 and A3		Adult NIBP Hose	ES-16-15	M1599A, M1599B, M1598B, 989803104341	YA25A06-10
		Neonate NIBP Hose	ESN-16-20	M1596B, M1597B	YA25A06-10N




Physio Control (See Medtronic)

Quick-Release NIBP Hose


Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Bionet: BMP300; Criticare: 506 series, 507EL/C, 507NJC, Datascope: Accutorr V, Duo, 2, Passport V, Spectrum, Trio, V12, V21; Edan; Goldway; Mindray T5; Nellcor; Ohmeda: CM5; Penlon; Philips: A1 and A3; Siemens; Spacelabs		Adult NIBP Hose	ES-15-15	200683-04-0003, 200683-04-0004, 509B-30-06259, 0683-00-0189-80, 705, 706, 708, 714-0060-00, 714- 714-0018-00	YA06A06-10
		Neonate NIBP Hose	ESN-15-20	708	YA06A06-10N

Siemens (See Drager and Quick-Release NIBP Hose)


Welch Allyn NIBP Hose (Also see Nihon Kohden)

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
MEK MP 700, MP 1000 and Welch Allyn Propaq Encore		Adult NIBP Hose	ES-09-09		YA09A09-10
		Neonate NIBP Hose	ESN-09-20		YA09A09-10N
Welch Allyn Atlas and Mennen monitors		Adult NIBP Hose	ES-04-09		YA04A04-10
		Neonate NIBP Hose	ESN-04-20		YA04A04-10N
Spot Vital Signs LXi		Adult NIBP Hose	ED-41-00	4500-32	YA41A41-10

Zoll NIBP Hose

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
M Series and E Series defibrillators		Adult NIBP Hose	ES-15-17	8000-0655, 8000-0662	YA15A17-10
		Neonate NIBP Hose	ESN-15-20		YA15A20-10

Tubing – no connectors

Models	Hose	Options	UNIMED/EVO	ORIGINAL OEM #	EVOTEC OEM
Tubing – no connectors 2.5m length		Single tubing 6mm - grey	SL0008A		SL0008A
		Double tubing 6mm - grey	SL0008B		SL0008B


Reusable NIBP Cuffs - **SOFT**

Description	Hose	Options	EVOTEC
Universal, reusable cuffs, sold individually. Connectors sold separately.		Adult (25 - 35cm) Single Hose	E1880S
		Adult (25 - 35cm) Double Hose	E1880D
		Large Adult (33-47cm) Single Hose	E1869S
		Large Adult (33-47cm) Double Hose	E1869D
		Small Adult (17-26cm) Single Hose	E1885S
		Small Adult (17-26cm) Double Hose	E1885D
		Pediatric (18-26cm) Single Hose	E1881S
		Pediatric (18-26cm) Double Hose	E1881D
		Infant (10-19cm) Single Hose	E1882S
		Infant (10-19cm) Double Hose	E1882D
		Neonate (6-11cm) Single Hose	E1883S
		Neonate (6-11cm) Double Hose	E1883D
		Thigh (46-66cm) Single Hose	E1884S
		Thigh (46-66cm) Double Hose	E1884D



Reusable NIBP Cuffs - **A-LOCK**

Description	Hose	Options	EVOTEC
Universal, reusable cuffs, sold individually. Connectors sold separately.		Adult (25 - 35cm) Single Hose	C850
		Adult (25 - 35cm) Double Hose	C840
		Large Adult (33-47cm) Single Hose	C850X
		Large Adult (33-47cm) Double Hose	C840X
		Pediatric (18-26cm) Single Hose	C850C
		Pediatric (18-26cm) Double Hose	C840C

Reusable NIBP Cuffs - **BLADDERLESS**

Description	Hose	Options	EVOTEC
Universal, reusable cuffs, sold individually. Connectors sold separately.		Adult (25 - 35cm) Single Hose	E1880NS
		Adult (25 - 35cm) Double Hose	E1880ND
		Large Adult (33-47cm) Single Hose	E1869NS
		Large Adult (33-47cm) Double Hose	E1869ND
		Pediatric (18-26cm) Single Hose	E1881NS
		Pediatric (18-26cm) Double Hose	E1881ND
		Infant (10-19cm) Single Hose	E1882NS
		Infant (10-19cm) Double Hose	E1882ND
		Neonate (6-11cm) Single Hose	E1883NS
		Neonate (6-11cm) Double Hose	E1883ND
Thigh (46-66cm) Single Hose	E1884NS		
Thigh (46-66cm) Double Hose	E1884ND		

Neonate Disposable NIBP Cuffs

Description	Hose	Options	EVOTEC
Disposable neonate NIBP cuffs in 5 sizes, sold in cases of 10 units.	 	Single Hose #1 (3-6cm)	E1781S-C05-10
		Double Hose #1 (3-6cm)	E1781D-C05-10
		Single Hose #2 (4-8cm)	E1782S-C05-10
		Double Hose #2 (4-8cm)	E1782D-C05-10
		Single Hose #3 (6-11cm)	E1783S-C05-10
		Double Hose #3 (6-11cm)	E1783D-C05-10
		Single Hose #4 (7-13cm)	E1784S-C05-10
		Double Hose #4 (7-13cm)	E1784D-C05-10
		Single Hose #5 (8-15cm)	E1785S-C05-10
		Double Hose #5 (8-15cm)	E1785D-C05-10

NIBP Connectors

Connector	Image Src	Description	UNIMED/EVO
A01 NIBP Connector		NIBP Replacement connector for various cuffs and hoses Replaces all expensive connectors between hose and cuff 1 to 1 connector	A01
A02 NIBP Connector		NIBP Female connector – connects with A05	A02
A03 NIBP Connector		NIBP luer-lock connector for use with Mennen (Old), Shchiller, Spacelabs (Old). Connects to A06	A03
A04 NIBP Connector		NIBP luer-style connector for use with Mennen (Old), Shchiller, Spacelabs (Old). Connects to A07.	A04
A05 NIBP Connector		NIBP connector for use with CAS Med, Drager, Kontron, Physio Control (Old), Schiller. Connects to A20. Generally used for neonate cuffs.	A05
A06 NIBP Connector		NIBP connector for use with BCI (Older versions), Datascope (Older versions). Connects to A03.	A06
A07 NIBP Connector		NIBP connector for use with BCI (Older versions), Datascope (Older versions). Connects to A04.	A07
A08 NIBP Connector		NIBP connector for use with Critikon, Invivo (requires 2), Kontron (requires 2), Mennen (requires 2), Nihon Kohden (requires 2), Welch Allyn. Connects to A09.and A09-P	A08
A09 NIBP Connector		NIBP threaded connector for use with Critikon, Invivo (requires 2), Kontron (requires 2), Mennen (requires 2), Nihon Kohden (requires 2), Welch Allyn. Connects to A08 .	A09

NIBP Connectors (Continued)

Connector	Image Src	Description	UNIMED/EVO
A09-P NIBP Connector		NIBP threaded connector for use with Critikon, Invivo (requires 2), Kontron (requires 2), Mennen (requires 2), Nihon Kohden (requires 2), Welch Allyn. Connects to A08 .	A09-P
A12 NIBP Connector		NIBP connector for use with Criticare, Datascope, Drager, Siemens, Physio Control, Mindray, Philips, Spacelabs. Connects to A15.	A12
A12-P NIBP Connector		NIBP connector for use with Criticare, Datascope, Drager, Siemens, Physio Control, Mindray, Philips, Spacelabs. Connects to A15.	A12-P
A15 NIBP Connector		NIBP bayonet-style connector for use with Schiller and Infinium One. Connects to A12.	A15
A16 NIBP Connector		NIBP connector for use with Philips. Connects to Philips monitors.	A16
A17 NIBP Connector		NIBP quick-connect connector for use with Datex Ohmeda, Marquette, Zoll M Series. Connects to A18.	A17
A18 NIBP Connector		NIBP quick-connect connector for use with Datex Ohmeda, Marquette, Zoll M Series. Connects to A17.	A18
A20 NIBP Connector		NIBP connector for use with neonate NIBP Hoses. Connects to A05.	A20
A21 NIBP Connector		NIBP connector for use with Pace Tech monitors.	A21

NIBP Connectors (Continued)

Connector	Connector	Connector	Image Src	UNIMED/EVO
A21S NIBP Socket		AS21 monitor socket. Connects to A21.		AS21S
A22 NIBP Connector		A22 NIBP connector for use with GE Datex AS/3, CS/3, S/5, Cardiocap 5, Cardiocap II, Eagle 4000, Solar, Tram monitors.		A22
A23 NIBP Connector		A23 NIBP connector for use with Marquette		A23
A24 NIBP Connector		A24 NIBP connector for use with IRIS and various other Connects to A24 S		A24
A24S NIBP socket		A24S NIBP connector socket for use with A24		A24S
A25 NIBP Connector		A25 NIBP connector for use with Dinamap monitors.		A25
A26S NIBP Socket		A26S NIBP Socket for use on monitor side Male Bayonet Socket to 5/32 in ID Barbed, Metal, L=29		A26S
A27S NIBP Socket		A27S NIBP Socket for use on monitor side Male Bayonet Socket to 5/32 in ID Barbed, Metal, L=36		A27S
A28 NIBP Connector		NIBP connector for use with Nihon Kohden monitors. NEW		A28

NIBP Connectors (Continued)

A31 NIBP Connector		NIBP connector for use with Colin monitors.	A31
A32 NIBP Connector		NIBP connector for use with Colin monitors.	A32
C875 Bulb		NIBP Bulb with valve for aneroids	C875
C888 Coiled tubing		NIBP Coiled tubing 8 foot for aneroids	C888