

EvoCare[®]



Temperature

Skin Probes

Reusable, 2.2m cable



Rectal/Esophageal

Reusable, 2.2m cable



Dual Probes

Double temp, 2.2m cable



Adapter Cables

For reusable & disposable probes, 2m cable



Disposable Temperature Probes



For adult, pediatric and infant warmers




Compatible,
reusable
probes

Adecon (See YSI 400)




Air Shields Reusable Temperature Probes (Infant warmers and Incubators)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
C-100		Pediatric Skin	EAS-15AS		EW0010A
		Pediatric Rectal Sensor	EAS-15AG		EAS-15AG
VERSALET C2000		Pediatric skin 1.2m			EW0045C







Air Shields Adaptor cable

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
VERSALET C2000 Fit Drager 4S and 4L temp sensors	MONITOR Connector 	Adaptor cable 0.67m			EW0045K

Air Shields Disposable sensors


Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
VERSALET C2000	MONITOR Connector 	Pediatric skin 1.2m		MU12520	EW0045K
Air Shields C2000 Isolette Incubators and Versalet 7700 Care Centers.		Premie Skin probe 5-l Box of 25	ETAS30-D-10 BOX OF 10	Replaces Air-Shields #83-104-15 Drager MU 12551 MU 12551	ET-10415 BOX OF 25
All Ohmeda Giraffe Incubators and the new Panda Warmer		Premie skin probe	ETOH30-D-10 BOX OF 10	Replaces Ohmeda 6600-0873-700 and 6600-0874-700 MU12520	ET-08730 BOX OF 25

Alaris Reusable Temperature Probes


Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
VITAL CHECK 4400/4410/4500 		Oral 4m cable	EVTOP-AR-001	3880L	EVTOP-AR-001
		Rectal 4m cable	EVTRP-AR-001	3882L	EVTRP-AR-001
Turbo Temp 2180C/2180CX01EE 		Oral 4m cable	EVTOP-AR-002	2885C	EVTOP-AR-002
		Rectal 4m cable	EVTRP-AR-002	2886C	EVTRP-AR-002
Turbo Temp 2185B/2185BX01EE 		Oral 4m cable	EVTOP-GE-001	2887B	EVTOP-GE-001
		Rectal 4m cable	EVTRP-GE-001	2888B	EVTRP-GE-001

AME (See Servocrib)



ATOM Reusable sensors (*Infant Warmers and Incubators*)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Infant Warmer Model V-3200 / V-505 Infant Incubator V2100G / V-2200 Transport Incubator V-707/ V-808		Pediatric Skin	ETAT-15PS		ETAT-15PS


Biolight Reusable sensors

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
V6 Vital Signs Anyview A2/ A3/ A5/ A8 		Adult Skin	ETBLT-AS		ETAT-15PS
		Pediatric Skin	ETBLT-PS		ETBLT-PS
		Adult general	ETBLT-AG		ETBLT-AG
		Pediatric General	ETBLT-PG		ETBLT-PG



Biolight Adaptor cable

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
V6 Vital Signs Anyview A2/ A3/ A5/ A8	Monitor Side 	0.3m	ETBLT-AD		ETBLT-AD
	Sensor Side 	For use with YSI400 Reusable probes			



Biolight Disposable Sensors

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
V6 Vital Signs Anyview A2/ A3/ A5/ A8		Skin 3m	ETBLT30-D10		ETBLT30-D10
		General 3m	ETBLT30-DAG-10		ETBLT30-DAG-10


Beijing BNT Reusable sensors (*Infant Warmers*)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Radiant Warmer BNT1000, BNT 2000, BNT 2000B		Pediatric Skin 3.0m	ET20-PS		ET20-PS
Radiant Warmer BNT1000, BNT 2000, BNT 2000B		Pediatric Skin 3.0m	ETBNT-15PS		ETBNT-15PS

Beijing BNT Disposable sensors (Infant Warmers)


Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Infant Radiant Warmer BNT 1000, 2000, 2000B		3m Direct connect Box of 10	ET2030-D-10		ET2030-D-10
Infant Radiant Warmer BNT 1000, 2000, 2000B		3m Direct connect Box of 10	ETBNT30-D-10		ETBNT30-D-10

Creative



Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
PC1000, 5000, 9000 UP6000		Skin 3m			EW0005A
		General 3m			EW0005B

Criticare (See YSI)

Datascope Reusable Temperature Probes (Also see YSI)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Datascope Passport V, V12, V21		Adult Skin Sensor	ETDT-AS	040-000057-00	ETDT-AS
		Adult Rectal Sensor	ETDT-AG	040-000055-00	ETDT-AG

Datascope Temperature Adapter Cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Datascope Passport V, V12, V21	To Monitor:  To YSI 400 Probe: 	0.3m Cable	ETDT-AD		ETDT-AD

Datascope Disposable Temperature Probes *(Also see YSI)*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Datascope Passport V, V12, V21		Adult Skin Sensor Box of 10	ETDT30-DAG-10		ETDT30-DAG-10
		Adult Rectal Sensor Box of 10	ETOH30-D-10		ETOH30-D-10

Datex *(See YSI)*

Dreager / Siemens Reusable Temperature Probes *(Also see YSI)*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Siemens Dreager Infinity Vista, Delta, Gamma, Kappa		Adult Skin Sensor	ETSM-AS	5204669	EW0008A
		Adult Rectal Sensor	ETSM-AG	5204644, 5204651	EW0008B

Dreager Temperature Adapters

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with disposable temperature probes	<p>To Monitor:</p> <p>To Philips Mono Probe:</p> <p>To GE Dual Pin Probe:</p>	Philips Mono Probe 0.3M	ETSM-AD	Compatible with Reusable YSI 400 series Philips probes	ETSM-AD
		GE Dual Pin Probe		11-CA-407-04, 11-CA-407-12	EW0008G

Dreager Disposable sensors

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
YSI 400 SERIES		Premie skin probe		Disposable YSI 400 series premie skin probe	EW0001E BOX OF 24
YSI 400 SERIES		ADULT Esophageal / rectal probe		Disposable YSI400 series Adult Esophageal / Rectal probe	EW0001f BOX OF 24

Drager Disposable sensors

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
YSI 400 SERIES		Premie skin probe		Disposable YSI 400 series premie skin probe With 3.5mm mono connector	EW0003E BOX OF 24
YSI 400 SERIES		ADULT Esophageal / rectal probe		Disposable YSI400 series Adult Esophageal / Rectal Probe With 3.5mm mono connector	EW0003F BOX OF 24

Drager *Radiant warmer / Incubator*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Infant Radiant warmer and Incubators		Premie skin probe Yellow	ETDG-D-20 BOX OF 20	Replaces Drager MX 11000	ET-11000 BOX OF 20
Infant Radiant warmer and Incubators		Premie skin probe		Replaces Drager MX 11000	EW0020E BOX OF 12




EDAN Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
EDAN M9		Adult Skin Sensor			EW0001A
		Adult Rectal Sensor			EW0001B







FANEM Reusable Temperature Probes (*Warmers and Incubators*)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Fanem Warmer and Incubator		Adult Skin Sensor			EW0066A
		Adult Rectal Sensor			EW0066B

Fisher & Paykel Reusable Temperature Probes (Humidifiers & anaesthetic)

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Air Temp probe for MR850 Humidifier		Air Temp		900MR869	900MR869
Air Temp probe for MR850 and HC550 Humidifier		Air Temp		900MR561	900MR561
MR850 Humidifier skin temp probe		Skin temp	ET700-AS		ET700-AS

GE Critikon Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
DINAMAP Pro200 / 400 V1 		Oral 4m cable	ETOP-AR-001	088012	ETOP-AR-001
		Rectal 4m cable	ETRP-AR-001	088016	ETRP-AR-001
Pro200/400 V2 ProCare200/400 		Oral 4m cable	ETOP-AR-002	2008774-001	ETOP-AR-002
		Rectal 4m cable	ETRP-AR-002	2008775-001	ETRP-AR-002
9300 Series DINAMAP XL 		Oral 4m cable	ETOP-GE-001	008975,008978	ETOP-GE-001
		Rectal 4m cable	ETRP-GE-001	008976	ETRP-GE-001

GE Marquette Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Carescape, Dash, Eagle, Tram, Solar, Mac- Lab, Cath-Lab		Adult Skin Sensor 3m	ETMQ-AS		ETMQ-AS
		Adult Rectal Sensor 3m	ETMQ-AG		ETMQ-AG
		Adult Skin Double Sensor 3m	ETMQ2-AS		ETMQ2-AS
7000/7010, GE Eagle, Solar, Tram, Tramscope, MAC-Lab Systems		Adult Skin Sensor 3m			EW0001A
		Adult Rectal Sensor 3m			EW0001B
		Pediatric Skin 3m			EW0001C
		Pediatric Rectal 3m			EW0001D

GE Marquette Temperature Adapters

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	
For use with disposable temperature probes	To Monitor: To Philips Mono Probe: To GE Dual Pin Probe: 	Philips Mono Probe	ETMQ-30-PH	2021700-001	ETMQ-30-PH
		GE Dual Pin Probe	ETMQ-30-AD	2021701-001	ETMQ-30-AD
For use with reusable YSI 400 temperature probes. 1 foot cable	To Monitor: To YSI 400 Probe: To DUAL 400 Probe: 	Single Temp	ETMQ-AD	2021700-001	ETMQ-AD
		Dual Temp YSI 400	ETMQ400-AD2	2016998-001	EW0015A
		Dual Temp YSI 700	ETMQ700-AD2	402015-004	EW0015A

**GE Marquette,
Infinium, MEDIX,
Mindray**

GE Marquette Disposable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with temperature adapter cables		Adult/Pediatric Skin, Case of 20	ETMQ-D-N	8001642	EW0001E
		Adult/Pediatric Esophageal/Rectal, Case of 20	ETMQ-DAG-N	8001643, 8001644, V3417, ER400-9	EW0001F
YSI 400 SERIES		Premie skin probe		Disposable YSI 400 series premie skin probe With 3.5mm mono connector	EW0003E BOX OF 24
YSI 400 SERIES		ADULT Esophageal / rectal probe		Disposable YSI400 series Adult Esophageal / Rectal Probe With 3.5mm mono connector	EW0003F BOX OF 24

Infinium *(See YSI)*

MEDIX Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
		Adult Skin Sensor		DIR-10	EW0010M



Mindray Reusable Temperature Probes *(Also see YSI)*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Mindray T5, T8 and iPM-9800		Adult Skin Sensor	ETMR-AS	0011-30-37393	EW0013A
		Adult Rectal Sensor	ETMR-AG	0011-30-37382	EW0013B


Mindray Temperature Adapters *(For use with disposable probes)*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with disposable temperature probes	To Monitor: To Philips Mono Probe: To GE Dual Pin Probe: 	Philips Mono Probe	ETMR-30-PH	011-30-37391, 0011-30-90444	ETMR-30-PH
		GE Dual Pin Probe	ETMR-30-AD		EW0013G

Mindray Temperature Adapters

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with reusable YSI 400 temperature probes	To Monitor:  To YSI 400 Probe: 	0.3M Cable	ETMR-AD		EW00013J

Mindray SMART TEMP



Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Accutor PLUS		Oral	ETOP-MR-001		EVTOP-MR-001
		Rectal	ETRP-MR-001		EVTRP-MR-001

Nihon Kohden *(See YSI)*

Ningbo David *Incubators and warmers*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
					EV-NB-2P
TI-2000		Skin sensor 3m			EV-NB-4P
HKN-90/HKN-93HKN-93A HKN-93B/HKN-9010/HKN-2000		Skin sensor 3m			EV-NB-6P

Ohmeda *Incubator*

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Incubators Giraffe Omnibed / Panda		Neonate Skin Sensor	ETOH-PS		EW-OH-S3P
		3meter			
NEW Panda warmer		Neonate Skin Sensor			EW-OH-SM
		3meter			

Ohmeda Disposable sensors

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
All Ohmeda Giraffe Incubators and the new Panda Warmer		Premie skin probe	ETOH30-D-10 BOX OF 10	Replaces Ohmeda 6600-0873-700 and 6600-0874-700 MU12520	ET-08730 BOX OF 25

OSEN Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
3000F		Adult Skin Sensor			EO3000F-S
		Adult Rectal Sensor			EO3000F-R


Philips Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with M1029A temperature module		Adult Skin Sensor	ETHP-AS	21078A	EW0003A
		Adult Rectal Sensor	ETHP-AG	21075A	EW0003B
For use with M1029A temperature module		Pediatric Skin Sensor	ETHP-PS		EW0003C
		Pediatric Rectal Sensor	ETHP-PG		EW0003D


Philips Temperature Adapters

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with disposable temperature probes	To Monitor: To Philips Mono Probe: To GE Dual Pin Probe:	Philips Mono Probe	ETHP-30-PH	21082A, 21082B, 98980310921	EW0003H
		GE Dual Pin Probe	ETHP-30-AD		EW0003G
For use with reusable YSI 400 temperature probes	To Monitor: To YSI 400 Probe:	0.3M CABLE	ETHP-AD		EW0003J


Philips Disposable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
For use with temperature adapter cables		Adult/Pediatric Skin, Case of 20	DHP-D-N0	21091A, 989803100951	EW0003E
		Adult/Pediatric Esophageal/Rectal, Case of 20	DHP-DAG-N0	21090A, 989803100941, M1837A, 989803105321	EW0003F



Servocrib Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Servo 2000 neonatal warmers		Adult/Pediatric Skin			EV-SERVO-S
		Adult/Pediatric Esophageal/Rectal, Case of 20			EV-SERVO-R



Siemens Reusable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Servo 2000 neonatal warmers		Adult / pediatric Skin	ETSM-AS	5204669	EW0008A
		Adult/ pediatric esophageal rectal	ETSM-AG	5204644	EW0008B

Siemens Adaptor cables

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Servo 2000 neonatal warmers		Adaptor cable to 400 series			EW0008G
			ETSM-AD		ETSM-AD


Siemens Disposable Temperature Probes

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
YSI 400 Series Disposable Sensors - use with EW0008G		Box of 24 Premie sensors with reflector			EW0001E
YSI 400 Series Disposable Sensors - use with EW0008G		Box of 24 Premie sensors			EW0001E


Sino Hero (See BLT)

Welch Allyn (See YSI)

YSI Reusable Temperature Probes - 400 SERIES

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
CASMed: 750; Criticare: 506, 507, 507E, 507ELC, 507S, 508, 1100, 2200, 8100 series; Datascope: Passport 2, Spectrum, Spectrum OR, Passport XG, Trio; Datex Ohmeda: S/3, CS/3, S/5, Cardiocap 1, Cardiocap 2, Cardiocap 5, Light Monitor, Satellite Plus; Drager: Delta, Vista, Gamma, appa; Invivo: MDE Prism; Mindray: PM 7000, PM 8000, PM 9000 Express, PM 9000, VS-800; Nihon Kohden: BSM-2300 Life Scope I, BSM-4100 Life Scope P, BSM-5100 Life Scope A, BSM-9510 Life Scope M, BSM-9800 Life Scope S.; Welch Allyn: Atlas		Adult Skin Sensor	ET2252-AS	409B, 401, 709B, 701, 01-02-0326, YSI-409A, 422, 5201335, M35373,	EW0001A
		Adult Rectal Sensor	ET2252-AG	420, 01-02-0327, 01-02-0379, 5201319, YSI-401, M1024247, 0206-02-0001	EW0001B
		Pediatric Skin Sensor	ET252-PS		EW0001C
		Pediatric Rectal Sensor	ET2252-PG		EW0001D

YSI Reusable Temperature Probes - 10K SERIES

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
		Adult Skin Sensor			EW0004A
		Adult Rectal Sensor			EW0004B
		Pediatric Skin Sensor			EW0004C
		Pediatric Rectal Sensor			EW0004D

YSI Reusable Temperature Probes -20K SERIES

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
		Adult Skin Sensor			EW0006A
		Adult Rectal Sensor			EW0006B


YSI Temperature Adapters for DISPOSABLE SENSORS

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
All monitors requiring YSI temperature sensors Disposable	<p>To Monitor:</p> <p>To Philips Mono Probe:</p> <p>To GE Dual Pin Probe:</p>	Jack Plug Philips	ETYSI-30-PH		EW0001H
		2 PIN ADAPTOR CABLE	ETYSI-30-AD	YSI 400/ 10K / 20K	EW0001G


YSI - DISPOSABLE SENSORS

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
All monitors requiring YSI temperature sensors Disposable		Philips Disposable skin Box of 20			EW0003E
		Philips disposable esophageal / rectal Box of 20			EW0003F


YSI - DISPOSABLE SENSORS

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
All monitors requiring YSI temperature sensors Disposable		YSI 400 series SKIN			EW0001E
		YSI 400 series Esophageal / Rectal 2			EW0001F


YSI – AIR SENSOR

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Continuous temperature measurement of gasses			ET100-GP		EW0001L



FOLEY Catheter - YSI 400 series

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Foley catheters with 400 series temperature sensor Silicone waterproof LATEX-FREE construction		14FR	ETMQ-CA-14F		ETF0014F
		16F	ETMQ-CA-16F		ETF0016F
		18F	ETMQ-CA-18F		ETF0018F

Tympanic temperature probe - Disposable - YSI 400 series

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Tympanic ear canal temperature measurement Dual foam push tube insertion		Adult	ETMQ-AE		ETT001A
		Pediatric	ETMQ-PE		ETT001P

Esophageal Stethoscope - Disposbale - YSI 400 series

Models	Connector	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Esophageal stethoscope with 400 series temperature sensor		9FR Pediatric	ETMQ-OP-9F		ETE0009F
		12F Adult	ETMQ-OP-12F		ETE0012F
		18FR Adult	ETMQ-OP-18F		ETE0018F
		24FR Adult	ETMQ-OP-24F		ETE0014F
Earpiece for Esophageal Stethoscope		Ear Piece for Stethoscope	ETSTP-EP		ET-EAR